

When extracting teeth, it can be difficult to position a straight dental elevator due to restricted access in some areas in the mouth.

We recommend using "winged" elevators that are either "Inside-Bent" or "Back-Bent" for those areas.

The image below shows elevating on the palatal aspect of the mesiopalatal crown-root segment of 108 in a 4 year old male neutered Papillon with advanced periodontal disease. A size 2 "Inside-Bent" winged elevator is being used.

To download the catalogue with a selection of "Inside-Bent" and "Back-Bent" winged dental elevators, go to the website www.cislak.com and click on Section 2 (Part 1).

